

THE EXTRACTION OF COPPER FROM COPPER CARBONATE (STEP II) - THE ELECTROLYSIS OF COPPER (II) SULFATE SOLUTION

Apparatus

Comboplate®
LED + crocodile clips
9V battery
Carbon electrodes (graphite pencil leads)




Eye
protection
must be worn
in this area

Method

1. Set up the apparatus as shown in the diagram below.
2. Leave the experiment to run for a few minutes and observe what takes place at both electrodes.

Chemicals

Copper(II) sulfate solution from Step I and



Results

| Electrode | Observations |
|-------------|--------------|
| Cathode (-) | |
| Anode (+) | |



**THE EXTRACTION OF
COPPER FROM COPPER
CARBONATE (STEP II) - THE
ELECTROLYSIS OF COPPER
(II) SULFATE SOLUTION**

Conclusions

Explain the observations made at both electrodes. _____

Give the equation for the production of copper at the negative electrode. _____

Explain why the Cu^{2+} ion is reduced in this reaction. _____

Explain why this reaction could not take place at the anode (+). _____

